



replacement parts industries, inc.

"The Alternate Source"

CALL 800-221-9723 or 818-882-8611 • FAX 818-882-7028

E-MAIL: techsupport@rpiparts.com • WEBSITE: www.rpiparts.com

ATV661 DRAIN VALVE

IMPORTANT INFORMATION REGARDING THE KNOB! READ BEFORE USING

The round-style knob on the RPI Drain Valve (RPI Part #ATV661) rotates clockwise to close the valve, while the OEM lever-style handle rotates counter-clockwise to close the valve.

- If installing the enclosed RPI drain valve with the Supplied round-style knob, it is important to replace the existing OEM drain rotation decal with the enclosed RPI decal.
- If installing the enclosed RPI drain valve and reusing the OEM lever-style handle, it is *not* necessary to replace the existing OEM drain rotation decal.