



# Castle/Getinge Bulk Sterilizers

Models: 3522, 3523, 3525, 3533, 3622, 3623 & 3633

## PM PACKAGES & SPARE PARTS PACKAGES

RPI PART #	OEM PART #	DESCRIPTION	(QTY) & RPI PART #'S INCLUDED IN KIT
CSP036	61301602085	<b>PM PKG (SOLENOID VALVES)</b> Models: 3522, 3523 & 3533	• Solenoid Valve Repair Kits: (1) AMK174, (2) AMK177, (1) AMK183 & (2) CSK048
CSP037	61301602090	<b>PM PKG (SOLENOID VALVES)</b> Models: 3525	• Diaphragm Replacement: (1) CSD051 • Solenoid Valve Repair Kits: (1) AMK174, (2) AMK177, (1) AMK183, (1) CSK048 & (1) CSK050
CSP038	61301602091	<b>PM PKG (FILTER &amp; SCREEN)</b> Models: 3522, 3523, 3525 & 3533	• Filter: (1) CSF052 • Screen Filter: (2) CSS053
CSP039	61301602087	<b>PM PKG (STM TRAP &amp; CHK VLV)</b> Models: 3525	• Check Valve Rebuild Kit: (3) CSK066 • Steam Trap: (2) CST059
CSP040	61301602092	<b>PM PKG (STM TRAP &amp; CHK VLV)</b> Models: 3525	• Check Valves: (1) CSV056, (1) CSV057 & (1) CSV058 • Check Valve Rebuild Kit: (1) CSK066
CSP068	61301602079	<b>PM PKG (STM TRAP &amp; CHK VLV)</b> Models: 3622 & 3633	• Check Valve Rebuild Kits: (1) CSK065 & (3) CSK066 • Steam Trap: (1) CST059 • Steam Trap Rebuild Kit: (1) CSK060
CSP069	No OEM Part # Available	<b>PM PKG (STM TRAP &amp; CHK VLV)</b> Models: 3522, 3523, 3525 & 3533	• Check Valve Rebuild Kit: (3) CSK066 • Steam Trap Rebuild Kit: (1) CSK067
CSP070	No OEM Part # Available	<b>PM PKG (STM TRAP &amp; CHK VLV)</b> Models: 3622 & 3633	• Check Valve Rebuild Kits: (1) CSK065 & (3) CSK066 • Steam Trap Rebuild Kit: (1) CSK060 & (1) CSK067
CSP041	61301602510	<b>SPARE PARTS PACKAGE</b> Models: 3525	• Check Valves: (1) CSV056, (1) CSV057 & (1) CSV058 • Check Valve Rebuild Kit: (3) CSK066 • Diaphragm Replacement: (1) CSD051 • Door Gasket: (1) CSG002 • Filter: (1) CSF052 • Solenoid Valve Repair Kits: (1) AMK174, (2) AMK177, (1) AMK183, (3) CSK048 & (1) CSK050 • Screen, Detergent: (1) CSS054 • Solid State Relay: (1) CSR028 & (1) CSR029
CSP042	61301602498	<b>SPARE PARTS PACKAGE</b> Models: 3522	• Check Valve Rebuild Kit: (3) CSK066 • Door Gasket: (1) CSG002 • Solenoid Valve Repair Kits: (1) AMK174, (2) AMK177 (1) AMK183 & (2) CSK048 • Solid State Relay: (1) CSR028 & (1) CSR029 • Steam Trap: (2) CST059
CSP043	61301602499	<b>SPARE PARTS PACKAGE</b> Models: 3523	• Check Valve Rebuild Kit: (3) CSK066 • Door Gasket: (1) CSG002 • Filter: (1) CSF052 • Solenoid Valve Repair Kits: (1) AMK174, (2) AMK177, (1) AMK183 & (3) CSK048 • Solid State Relay: (1) CSR028 & (1) CSR029 • Steam Trap: (3) CST059
CSP044	61301602500	<b>SPARE PARTS PACKAGE</b> Models: 3533	• Check Valve Rebuild Kit: (3) CSK066 • Door Gasket: (1) CSG002 • Filter: (1) CSF052 • Solenoid Valve Repair Kits: (1) AMK174, (2) AMK177, (1) AMK183, (1) CSK018 & (2) CSK048 • Solid State Relay: (1) CSR028 & (1) CSR029 • Steam Trap: (2) CST059
CSP045	61301602501	<b>SPARE PARTS PACKAGE</b> Models: 3522	• Check Valve Rebuild Kit: (3) CSK066 • Door Gasket: (1) CSG001 • Solenoid Valve Repair Kits: (1) AMK174, (2) AMK177, (1) AMK183 & (2) CSK048 • Solid State Relay: (1) CSR028 & (1) CSR029 • Steam Trap: (2) CST059
CSP046	61301602502	<b>SPARE PARTS PACKAGE</b> Models: 3523	• Check Valve Rebuild Kit: (3) CSK066 • Door Gasket: (1) CSG001 • Filter: (1) CSF052 • Solenoid Valve Repair Kits: (1) AMK174, (2) AMK177, (1) AMK183 & (3) CSK048 • Solid State Relay: (1) CSR028 & (1) CSR029 • Steam Trap: (3) CST059
CSP047	61301602503	<b>SPARE PARTS PACKAGE</b> Models: 3533	• Check Valve Rebuild Kit: (3) CSK066 • Door Gasket: (1) CSG001 • Filter: (1) CSF052 • Solenoid Valve Repair Kits: (1) AMK174, (2) AMK177, (1) AMK183, (1) CSK018 & (2) CSK048 • Solid State Relay: (1) CSR028 & (1) CSR029 • Steam Trap: (2) CST059