




## MORE NEW PARTS TO FIT GENDEX GXP™ FILM PROCESSORS

**Ready-Made Clear Tubing Kits ... You Can See the Solution Flow Inside the Tubing!**




RPI Part #GXK026  
OEM Part #112-0623G1

**INLET TUBING KIT (RED - 3/8" Tubing)**  
• Includes 6' of 3/8" ID PVC Clear Tubing (RPI Part #RPT646), (2) 3/8" OD Red Spiral Insert (RPI Part #GXS032) & (2) Kwik™ Clamp (RPI Part #RPC636)

RPI Part #GXK029  
OEM Part #112-0622G1

**OUTLET TUBING KIT (RED - 1/2" Tubing)**  
• Includes 6' of 1/2" ID PVC Clear Tubing (RPI Part #RPT647), (2) 1/2" OD Red Spiral Insert (RPI Part #GXS035) & (2) Kwik™ Clamp (RPI Part #RPC286)




RPI Part #GXK027  
OEM Part #112-0623G2

**INLET TUBING KIT (BLUE - 3/8" Tubing)**  
• Includes 6' of 3/8" ID PVC Clear Tubing (RPI Part #RPT646), (2) 3/8" OD Blue Spiral Insert (RPI Part #GXS033) & (2) Kwik™ Clamp (RPI Part #RPC636)

RPI Part #GXK030  
OEM Part #112-0622G2

**OUTLET TUBING KIT (BLUE - 1/2" Tubing)**  
• Includes 6' of 1/2" ID PVC Clear Tubing (RPI Part #RPT647), (2) 1/2" OD Blue Spiral Insert (RPI Part #GXS036) & (2) Kwik™ Clamp (RPI Part #RPC286)



RPI Part #GXK028  
OEM Part #112-0623G3

**INLET TUBING KIT (CLEAR- 3/8" Tubing)**  
• Includes 6' of 3/8" ID PVC Clear Tubing (RPI Part #RPT646), (2) 3/8" OD Black Spiral Insert (RPI Part #GXS034) & (2) Kwik™ Clamp (RPI Part #RPC636)

RPI Part #GXK031  
OEM Part #112-0622G3

**OUTLET TUBING KIT (CLEAR - 1/2" Tubing)**  
• Includes 6' of 1/2" ID PVC Clear Tubing (RPI Part #RPT647), (2) 1/2" OD Black Spiral Insert (RPI Part #GXS037) & (2) Kwik™ Clamp (RPI Part #RPC286)

### Create Your Own Clear Tubing Kits ... For Those Custom Jobs!

RPI Part #GXS032  
SPIRAL INSERT (3/8" RED)



SOLD IN 36" LENGTHS

RPI Part #GXS035  
SPIRAL INSERT (1/2" RED)


RPI Part #GXS033  
SPIRAL INSERT (3/8" BLUE)



SOLD IN 36" LENGTHS


RPI Part #GXS036  
SPIRAL INSERT (1/2" BLUE)

RPI Part #GXS034  
OEM Part #407-0156P1  
SPIRAL INSERT (3/8" BLACK)



SOLD IN 36" LENGTHS

RPI Part #GXS037  
OEM Part #407-0156P2  
SPIRAL INSERT (1/2" BLACK)



- Flexible clear PVC
- Sold by the foot (Reel not included)
- For Outet and Inlet Pump Connections

RPI Part #RPT646  
OEM Part # (No OEM Part # Available)  
PVC CLEAR TUBING (3/8" ID x 1/2" OD)

RPI Part #RPT647  
OEM Part # (No OEM Part # Available)  
PVC CLEAR TUBING (1/2" ID x 11/16" OD)

RPI Part #RPC286  
KWIK™ CLAMPS (6/pkg)  
• Fits 11/16" OD Tubing  
• White Nylon




RPI Part #RPC636  
KWIK™ CLAMPS (6/pkg)  
• Fits 1/2" OD Tubing  
• White Nylon



RPI Part #GXK040  
OEM Part #406-0015P1

**HINGE KIT**  
• 2 per package  
• Includes (2) mounting screws



Fits: Upper cover to lower base

RPI Part #RPS639  
OEM Part # (No OEM # Available)

**SILICONE GASKET MAKER 598**  
• 5 units/pkg; .05 fl. oz. per unit  
• Provides a sealant for pump tubing connections




RPI Part #RPC637  
OEM Part #407-0159P5

**KWIK™ CLAMPS (6/pkg)**  
Fits: 3/8" OD Tubing



RPI Part #RPT652  
OEM Part #273-0010P1

**TUBING (1/4" TYGON®)**  
• 1/4" ID x 3/8" OD; Black  
Fits: Pump and Pump Connections



DOES NOT INCLUDE REEL




replacement parts industries, inc.

CALL (800) 221-9723 • FAX (818) 882-7028 • EMAIL [order@riparts.com](mailto:order@riparts.com) • WEBSITE [www.riparts.com](http://www.riparts.com)

The above parts are manufactured by Replacement Parts Industries, Inc. to fit Gendex equipment. All product names used in this document are trademarks or registered trademarks of their respective holders. WEB Flyer 06/07

## ... Even More RPI Parts to fit Gendex GXP™ Film Processors!

- COOLING FAN
- DRIVE BELT
- HEATER FAN
- HINGE KIT
- LEVELING FOOT
- MAIN FUSE & HOLDER
- MOTOR BRUSH & BRUSH CAP
- O-RING (FITS STANDPIPE)
- PC BOARD FUSES (F1-4)
- PUSH BUTTON SWITCH (RED)
- PUSH BUTTON SWITCH (WHITE)
- SWITCH COVER



An ISO 9001 Certified Company

PRESORTED  
STANDARD  
U.S. POSTAGE  
PAID  
VAN NUYS, CA  
PERMIT #2399

• Address Service Requested

P.O. Box 5019 • Chatsworth, California • 91313-5019

replacement parts Industries, Inc.



## Celebrating More New Parts from RPI to Fit Gendex Film GXP™ Processors including Tubing!

### 3 Very Good Reasons Why RPI Inlet and Outlet Tubing Has the Advantage!

#### #1 Inlet & Outlet Tubing is CLEAR!

You can see solution flow inside the tubing!


#### #2 Color Coded Spiral Inserts – Red, Blue and Black!

Makes it easy to match up tubing with hose interface.


#### #3 All Parts in Kits are also Sold Separately!

Create tubing for those custom fit jobs.

3/8" ID Inlet Tubing – GXK026  
1/2" ID Outlet Tubing – GXK029



3/8" ID Inlet Tubing – GXK027  
1/2" ID Outlet Tubing – GXK030



3/8" ID Inlet Tubing – GXK028  
1/2" ID Outlet Tubing – GXK031

