





INTRODUCING NEW PARTS TO FIT GENDEX GXP™ FILM PROCESSORS




RPI Part #GXF012
OEM Part #652-0012P1
HEATER FAN (115V)
 • 4.69" square x 1.5" wide
Fits: Heater/Fan Assembly




RPI Part #GXF011
OEM Part #122-0096G1
COOLING FAN (115V)
 • 3.15" square x 1" wide
Fits: Chassis



RPI Part #GXB009
OEM Part #420-0007P1
DRIVE BELT
 • 190XL025 Power Grip Timing Belt
Fits: Developer, fixer, wash and dryer racks




RPI Part #GXS014
OEM Part #632-0045P1
PUSH BUTTON SWITCH (RED)
 • Endo switch
 • Includes Switch Cover (RPI Part #GXC015)
Fits: Switch panel




RPI Part #GXS013
OEM Part #632-0046P1
PUSH BUTTON SWITCH (WHITE)
 • Run/Stand-by and On/Off switches
 • Includes Switch Cover (RPI Part #GXC015)
Fits: Switch panel




RPI Part #GXC015
OEM Part #433-0060P1
SWITCH COVER
 • Transparent
Fits: All push button switches and lamp



RPI Part #GXB004
OEM Part #652-0015P1/Brush
MOTOR BRUSH (2/pkg)
 • Cap sold separately (RPI Part #RPC038)
Fits: Brush holder assembly on main drive motor



RPI Part #RPH638
OEM Part # 646-0047P1
FUSE HOLDER (MAIN)
 • 3/16" terminals
 • Fits standard 1/4" x 1-1/4" fuse
 • Includes mounting hardware
Fits: Rear panel



RPI Part #GXF024
OEM Part #112-0629G1
LEVELING FOOT
 • 4 per package
 • Non-skid PVC base
Fits: Bottom frame



RPI Part #RPC038
OEM Part #652-0015P1/Cap
BRUSH CAP (6/pkg)
 • Brush sold separately (RPI Part #GXB004)
Fits: Main drive motor




RPI Part #RPF070
OEM Part #646-0048P1
FUSE (MAIN)
 • 12A, 250V; 1/4" x 1-1/4"
 • Ceramic
Fits: Rear panel

PC BOARD FUSES: F1 TO F4

RPI Part #	OEM Part #	Position & Description
RPF045	640-0054P16	F1; 1A, 250V
RPF046	46-1700-21P34	F2; 1-1/2A, 250V
RPF053	46-1700-21P33	F3; 8A, 250V
RPF512	No OEM Part # Available	F4; 1A, 250V

TO HELP FIX THOSE ANNOYING LEAKS!

RPI Part #RPS639
OEM Part # (No OEM # Available)
SILICONE GASKET MAKER 598
 • 5 per package
 • .05 fl. oz. per unit
 • Provides a sealant for pump tubing connections



RPI Part #RPC637
OEM Part # 407-0159P1
KWIK™ CLAMP (6/pkg)
 • Range: 23/64"-25/64"
 • Use with pump tubing



RPI Part #RPO360
OEM Part #404-0036P2
O-RING (6/pkg)
 • .301" x .441" x .070"
 • Material: EPDM
Fits: Standpipe




LOOK FOR MORE PARTS TO FIT GENDEX FILM PROCESSORS AND X-RAY EQUIPMENT LATER THIS SUMMER!



RPI introduces a new line ...
parts to fit
GENDEX
GXP™ FILM PROCESSORS


...and many more to come later this summer!



**CUSTOMIZE YOUR WORK SPACE
 WITH RPI POWER CORDS!**

*Depending on where your outlet is located ...
 choose the power cord that best fits your needs!*

Right Angle Plug




RPI Part #RPC582
OEM Part #644-0009P1

INDUSTRIAL GRADE POWER CORD

- 14/3 SJT, 105°C; 15A@125V
- NEMA 6-15P; EN60320-C13R Right Angle

Straight Plug



RPI Part #RPC262


OEM Part # (No OEM # Available)

INDUSTRIAL GRADE POWER CORD

- 14/3 SJT, 60°C; 15A@125V
- NEMA 5-15P; EN60320-C13

**NEED THE
 POWER MODULE?**

*Check the power module.
 Is it time to replace it?*



RPI Part #HBM011
OEM Part #651-0068P1

POWER MODULE

- 15A@250V