



NEW! PCB's & SWITCHES TO FIT PELTON & CRANE CHAIRMAN, CHAIRMAN 5000 & COACHMAN DENTAL CHAIRS



CHAIRMAN


RPI Part #PCS720
New OEM Part #007899

EXIT SWITCH ASSEMBLY

- Normally Open
- 3/4A @ 125VAC
- Black button

Fits: Chair back assembly

Models: Chairman




RPI Part #PCS719
New OEM Part #007898

OPERATE SWITCH ASSEMBLY

- Normally Open
- 3/4A @ 125VAC
- Red button

Fits: Chair back assembly

Models: Chairman




RPI Part #PCS713
OEM Part #007431

ROLLER MICRO SWITCH

- 5A @ 250 VAC

Fits: Seat wiring assembly

Model: Chairman




RPI Part #PCS714
New OEM Part #007506

FORMED LEVER SWITCH

- 5A @ 120VAC

Fits: Base, lift & attachment assembly

Models: Chairman




RPI Part #PCS715
OEM Part #007596

TRAVERSE LIMIT SWITCH

- 5A @ 125/250VAC

Fits: Base, lift & attachment assembly

Models: Chairman




RPI Part #PCP721
OEM Part #007595

FOOT SWITCH PEDAL

- Black plastic

Fits: Chair base foot switch assembly

Models: Chairman




RPI Part #PCS716
OEM Part #007429

TRAVERSE SWITCH

- Momentary contact (On/Off/On)
- 6A @ 125 V; 3A @ 250 V
- Snap in mounting with low profile paddle (dark grey paddle)

Fits: Back assembly and foot switch assembly

Model: Chairman




RPI Part #PCS717
OEM Part #007428

TILT SWITCH

- Momentary contact (On/Off/On)
- 6A @ 125 V; 3A @ 250 V; Snap in mounting with low profile paddle (white paddle)

Fits: Back assembly (Chairman) and Foot Switch Assembly (Chairman 5090)

Model: Chairman & Chairman 5090




RPI Part #PCS718
OEM Part #007427

LIFT SWITCH

- Momentary contact (On/Off/On)
- 6A @ 125 V; 3A @ 250 V
- Snap in mounting with low profile paddle (black paddle)

Fits: Back assembly

Model: Chairman



COACHMAN & COACHMAN II


RPI Part #PCC728
OEM Part #94-416-43

LIFT CAPACITOR (56-60 µf)

- Plastic; Electrolytic medium

Fits: Base motor/lift motor

Models: Chairman 5000 & Coachman



RPI Part #PCC734
OEM Part #77-31-904/019518

CAPACITOR (21-25 µf)

- Plastic; Electrolytic medium

Fits: Incline/recline and upper traverse

Models: Chairman 5000 & Coachman




RPI Part #PCB723
New OEM Part #018326

SWITCH BACK PCB ASSEMBLY

Fits: Chair back assembly

Models: Coachman




RPI Part #PCS724
OEM Part #018420

MICRO SWITCH

- 5A @ 250 VAC

Fits: Chair back assembly

Model: Coachman




RPI Part #PCS722
OEM Part #023944

FOOT SWITCH MICROSWITCH

- 15A @ 125VAC

Fits: 4-way foot switch assembly

Models: Coachman & Coachman II




RPI Part #PCS725
OEM Part #018956

MOTOR LIMIT SWITCH

- 15A @ 120VAC

Fits: Base assembly

Models: Coachman & Coachman II



replacement parts industries, inc.

CALL (800) 221-9723 • FAX (818) 882-7028 • EMAIL order@rpiparts.com • WEBSITE www.rpiparts.com

The above parts are manufactured by Replacement Parts Industries, Inc. to fit Pelton & Crane equipment. All product names used in this document are trademarks or registered trademarks of their respective holders. WEBSITE 11/04



NEW! PCB's & SWITCHES TO FIT PELTON & CRANE CHAIRMAN, CHAIRMAN 5000 & COACHMAN DENTAL CHAIRS



CHAIRMAN 5000


RPI Part #PCC734
OEM Part #77-31-904/019518

CAPACITOR (21-25 µf)
• Plastic, Electrolytic medium
Fits: Incline/recline and upper traverse
Models: Chairman 5000 and Coachman




RPI Part #PCC728
OEM Part #94-416-43

LIFT CAPACITOR (56-60 µf)
• Plastic; Electrolytic medium
Fits: Base motor/lift motor
Models: Chairman 5000 & Coachman




RPI Part #PCB736
New OEM Part #94 34 424

UP INTERCONNECT (PCB)
Fits: Incline/recline controls
Models: Chairman 5000 & 5010




RPI Part #PCS735
OEM Part #77 40 947

A/C POWER SWITCH
• 16A @ 125/250VAC
• Use only with 240V models
Fits: Inner base assembly
Models: Chairman 5000 (240V models only)



RPI Part #PCS731
OEM Part #15-25-737


CORRECTION LIMIT SWITCH
• 10A @ 125 VAC
Fits: Lift motor
Model: Chairman 5000



CHAIRMAN 5010


RPI Part #PCB736
New OEM Part #94 34 424

UP INTERCONNECT (PCB)
Fits: Incline/recline controls
Models: Chairman 5000 & 5010



RPI Part #RPF367
OEM Part #18-81-163


FUSE (2A, 250V)
• 5 per package
• Size: 5 mm x 20 mm (.197" x .796")
• MDL Glass, Fast-Acting
Fits: Inner base assembly (240V models only)
Models: Chairman 5000 & 5010



CHAIRMAN 5090

RPI Part #PCB733
New OEM Part #30 06 629

LIMIT SWITCH (PCB)
Fits: Incline/recline and upper traverse
Models: Chairman 5090




RPI Part #PCA737
OEM Part #18 99 702

FOOT SWITCH & CABLE ASSEMBLY
• Momentary contact (On/Off/On)
• 6A @ 125 V; 3A @ 250 V; Snap in mounting
with low profile paddle (white paddle)
Fits: Inner base
Models: Chairman 5090



RPI Part #PCS717
OEM Part #007428

TILT SWITCH
• Momentary contact (On/Off/On)
• 6A @ 125 V; 3A @ 250 V; Snap in mounting
with low profile paddle (white paddle)
Fits: Back assembly (Chairman) and Foot Switch
Assembly (Chairman 5090)
Model: Chairman & Chairman 5090



replacement parts industries, inc.

CALL (800) 221-9723 • FAX (818) 882-7028 • EMAIL order@rpiparts.com • WEBSITE www.rpiparts.com

The above parts are manufactured by Replacement Parts Industries, Inc. to fit Pelton & Crane equipment. All product names used in this document are trademarks or registered trademarks of their respective holders. WEBSITE 11/04