



# NEW PARTS TO FIT MIDMARK•RITTER TABLES & CHAIRS MODELS EVOLUTION 75 and SERIES 100, 300 & 400



**RPI Part #MIH087**  
**OEM Part #002-0126-00**

**RETURN HOSE KIT**  
• Length: 24.25"; Neoprene  
**Fits:** Tilt to Back  
**Models:** See Chart Below

**RPI Part #MIH085**  
**OEM Part # 002-0124-00**


**RETURN HOSE KIT**  
• Length: 38.25"; Neoprene  
**Fits:** Foot to Tilt  
**Models:** See Chart Below

**RPI Part #MIH084**  
**OEM Part #002-0123-00**

**POWER HOSE KIT**  
• Length: 44.50"; Neoprene  
**Fits:** Foot to Tilt  
**Models:** See Chart Below

**RPI Part #MIC088**  
**OEM Part #016-0240-00**

**LOCKING GAS CYLINDER**  
**Fits:** Back rest assembly  
**Models:** See Chart Below



## RPI PARTS TO FIT MIDMARK•RITTER TABLES & CHAIRS MODELS EVOLUTION 75 and SERIES 100, 300 & 400

OEM PART #	RPI PART #	DESCRIPTION	Evolution 75	101	102	104	105	106	107	110	111	112	113	114	115	116	117	119	311	317	319	404	411	413	415	419
002-0001-00	MIC063	BASE CYLINDER KIT	1				24	25		1	27	1		1		1	1	1	1	1						
002-0002-00	MIC064	TILT CYLINDER KIT	1							1	27			1		1		1	1	1						
002-0003-00	MIC065	BACK/FOOT CYLINDER KIT	38				24			1	27			1		1	1	1	44		44					
002-0012-00	MIH068	RETURN HOSE KIT									3			3												
002-0014-00	MIH067	POWER HOSE KIT									28			3												
002-0021-00	MIH070	RETURN HOSE KIT									29															
002-0031-00	MIH069	HOSE ASSEMBLY KIT						4			3			3												
002-0032-00	MIH071	HOSE ASSEMBLY KIT						26			29															
002-0041-00	MIR017	TIME DELAY RELAY KIT	1					1		1	1	1	1	1		1	1	1	1	1						
002-0043-00	MIC020	CAPACITOR KIT						4		1	3			3												
002-0044-00	MIC019	CAPACITOR	39				12	15	1		13	1	1	13	1	1	1	1	1	1	1	1				
002-0085-00	MIS021	REED SWITCH ASSEMBLY													5											
002-0086-00	MIS022	REED SWITCH ASSEMBLY													1											
002-0087-00	MIS023	FOOT CONTROL SWITCH													1											
002-0101-00	MIS015	FOOT CONTROL SWITCH					7	1	1								1							1		
002-0113-00	MIS013	HEATER SWITCH (RED)		2		2	2																			
002-0114-00	MIS014	HEATER SWITCH (RED)		11	1	11	11		1										1							
002-0117-00	MIH001	PUMP HOSE KIT	1					19			30	1		20		1	1	1	1		1					
002-0118-00	MIH002	RETURN HOSE KIT	40								31	1		20		1	1	1	1							
002-0119-00	MIH003	POWER HOSE KIT	1								30	1		20		1	1	1		1	1					
002-0120-00	MIH004	RETURN HOSE KIT	1								30	1		20		1	1	1	1	1	1					
002-0121-00	MIH005	POWER HOSE KIT	1								30	1		20		1	1	1	1	1	1					
002-0122-00	MIH006	RETURN HOSE KIT	1								30	1		20		1	1	1	1	1	1					
002-0123-00	MIH084	RETURN HOSE KIT	1								30			20		1		1	1		1					
002-0124-00	MIH085	POWER HOSE KIT	1								32			20		1		1	1		1					
002-0125-00	MIH007	POWER HOSE KIT	1								30	1		20		1	1	1	1	1	1					
002-0126-00	MIH087	RETURN HOSE KIT	41								33			34			35	36	1	1	1	1				
002-0171-00	MIR024	CONTROL PANEL RELAY													16											
002-0178-00	MIH008	RETURN HOSE KIT	42								32			23		18	1	18								
014-0056-00	RPF384	HYDRAULIC FLUID					1	1	1	1	1	1		1	1	1	1	1								
015-0374-00	MIR018	AUTO RETURN & INTLCK RELAY									14	1	1	14							1					
015-0376-00	MIS012	STOP SWITCH (RED)									14	1	1	14							1					
015-0381-00	MIS009	RETURN LIMIT SWITCH									21	6	6	21				10								
015-0408-00	MIS025	ROCKER SWITCH																1			1					
015-0421-00	MIS016	RETURN LIMIT SWITCH									22	17	17	22				9			1					
015-0422-00	MIS026	RETURN LIMIT SWITCH																8			1					
015-0424-00	MIS011	AUTO RETURN SWITCH (GREEN)									14	1	1	14				1			1					
015-0437-00	DEC610	CAPACITOR	43																			1		1	1	
016-0240-00	MIC088	LOCKING GAS CYLINDER																				1				

**FOOTNOTES**

1 - Fits all Serial #'s.  
2 - Prior to Serial #14536.  
3 - Prior to Serial #37420.  
4 - Prior to Serial #37941.  
5 - Prior to Serial #A1500.  
6 - Prior to Serial #AF1167 & AG1046.  
7 - Serial #H1020 & G1325.  
8 - Tilt From Serial #AJ1277.  
9 - Base From Serial #AJ1273 & Tilt from Serial #AJ1273.  
10 - Base Prior to Serial #AJ1273 & Tilt prior to Serial #AJ1273.  
11 - From Serial #14536 and above.  
12 - From Serial #19053.  
13 - From Serial #37420.  
14 - From Serial #K3770 & L1228.  
15 - From Serial #37941.  
16 - From Serial #A1500.  
17 - From Serial #AF1167 & AG1046.  
18 - From Serial #AN1875 & AZ1041.  
19 - From Serial #J1050 thru present; CC1000 thru CC1568.  
20 - From Serial #K2726 & L1139.  
21 - From Serial #K3770 & L1228 to K5010 & L1342.  
22 - From Serial #K5010 & L1342.  
23 - From Serial #K8179 & L1859.  
24 - From Serial #G1000, H1000, V1000, CA1000, CF1000, & CG1000 thru present.  
25 - From Serial #37940 thru present; J1000 thru present; CC1000 thru present.  
26 - From Serial #37940 thru present; J1000 thru J1049.  
27 - Prior to Serial #37420 thru present; from BX1000 thru present; CA1000 thru present.  
28 - Prior to Serial #37420 thru K2725.  
29 - From Serial #37420 thru K2725.  
30 - From Serial #K2726, L1139, BX1000; CA1000 thru present.  
31 - From Serial #K2726 thru present; from BX1000 thru BX2310; CA1000 thru CA1506.  
32 - From Serial #K8179 thru present; BX1000 thru BX1739; CA1000 thru CA1220.  
33 - From Serial #K2726 thru K8178; from BX1740 and CA1221 thru present.  
34 - From Serial #L1139 thru present.  
35 - From Serial #FE1000 thru present.  
36 - From Serial #AN1000 thru AN1874; from BY2006 thru BY2624; from CA1191 thru CA1506; from DF2695 thru present.  
38 - Both foot and back cylinder for Serial #DF 1000 thru DF 2694; All other Serial #'s used as foot cylinder ONLY.  
39 - From Serial #DF 1000 thru present; from EM 1000 thru present; from EY 1000 thru present; from DZ 1000 thru present; from KK 1000 thru present; from KL 1000 thru present.  
40 - From Serial #DF1000 thru present; from EM1000 thru EM4519; from EN1000 thru EN1081.  
41 - From Serial #DF2695 thru present; from EM1000 thru present; from EN1000 thru present; from EY1000 thru present; from EZ1000 thru present; from DZ1000 thru present; from KK1000 thru present; from KL1000 thru present.  
42 - From Serial #DF1000 thru DF2694.  
43 - From Serial #EN1000 thru present; from EZ1000 thru present.  
44 - Foot Cylinder only.