

# RPI is Your Best Source for Table Top Sterilizer Parts!



- 28 different PM Kits for table top sterilizers serviced most
- Outstanding Technical Service – FREE via phone, fax and e-mail
- Troubleshooting Guides and PM Checklists – FREE downloads from [www.rpiparts.com](http://www.rpiparts.com)

## PM Kit to fit AMSCO Eagle® 10 & 10 Plus




RPI Part #AMK234  
OEM Part #  
P426637-261  
(Gasket Only)

### STERILIZER PM KIT

- Includes Door Gasket (RPI Part #AMG161), Water Guard (RPI Part #AMG172), Filter (Fill Line) (RPI Part #PCF009), and a PM Check List

Models: Eagle 10 and Eagle 10 plus



## PM Kits to fit Midmark® • Ritter® 7, M7, M9 & M11




RPI Part #MIK134  
OEM Part #  
002-0243-00  
(Gasket Only)

### STERILIZER PM KIT

- Includes Fill/Vent Mesh Chamber Filter (RPI Part #MIF062), Door Gasket (RPI Part #RCG019), and a PM Check List

Models: M7 and 7




## PM Kits to fit SCICAN STATIM® 2000 & 5000

RPI Part #SCK035  
OEM Part #  
01-100028S  
(Gasket Only)

### STERILIZER PM KIT

- Includes Cassette Seal (RPI Part #SCS001), Liquid Soap (RPI Part #RPS287), Compressor Filter (RPI Part #SCF002), Biological Filter (RPI Part #SCF032 for S/N 2101DK0274 and above), Water Filter (In-Line) (RPI Part #SCF034), and a PM Check List

Model: 2000




RPI Part #SCK036  
OEM Part #  
01-101649S  
(Gasket Only)

### STERILIZER PM KIT

- Includes Cassette Seal (RPI Part #SCS029), Liquid Soap (RPI Part #RPS287), Air Compressor Filter (RPI Part #SCF031 for use with MEDO Compressor only), Biological Filter (RPI Part #SCF032), Water Filter (In-Line) (RPI Part #SCF034), and a PM Check List

Model: 5000




RPI Part #MIK072  
OEM Part #  
002-0361-00  
(Door & Dam Gaskets only)

### STERILIZER PM KIT

- Includes Mesh Style Filters (RPI Part #MIF062), Sintered Bronze Filter (RPI Part #MIF038), Door Springs (RPI Part #MIS076), Door Gasket (RPI Part #MIG027), Dam Gasket (RPI Part #MIG035), and a PM Check List

Models: M9 & M9D




RPI Part #MIK080  
OEM Part #  
002-0504-00  
(Door & Dam Gaskets only)

### STERILIZER PM KIT

- Includes Mesh Style Filters (RPI Part #MIF062), Sintered Bronze Filter (RPI Part #MIF038), Door Springs (RPI Part #MIS076), Door Gasket (RPI Part #MIG028), Dam Gasket (RPI Part #MIG036), and a PM Check List

Models: M11 & M11D



replacement parts industries, inc.

CALL (800) 221-9723 • FAX (818) 882-7028 • EMAIL [order@rpiparts.com](mailto:order@rpiparts.com) • WEBSITE [www.rpiparts.com](http://www.rpiparts.com)

The above parts are manufactured by Replacement Parts Industries, Inc. to fit the equipment noted. All product names used in this document are trademarks or registered trademarks of their respective holders. Web 01/11

STERILIZER PM KITS – AMERICAN STERILIZER, MIDMARK® • RITTER®, PELTON & CRANE, SCICAN, & TUTTNAUER®

SECTION

10

PAGE


116

UPDATED

01/11



## PM Kits to fit Tuttnauer EZ9, EZ10, EZ10k, 1730, 2340, 2540 & 3870 series




**RPI Part #TUK121**  
**OEM Part #02610020**  
(Door Gasket Only)

### STERILIZER PM KIT

- Includes Fill/Vent Mesh Chamber Filter (RPI Part #MIF062), Door Bellows (RPI Part #TUB019), Rubber Boot For Door Switch (RPI Part #TUB072), Door Gasket (RPI Part #TUG001), and Air Jet Valve (Black Top) (RPI Part #TUJ034), and a PM Check List

**Models:** 1730E/M




**RPI Part #TUK126**  
**OEM Part #02610118**  
(Door Gasket Only)

### STERILIZER PM KIT

- Includes Fill/Vent Mesh Chamber Filter (RPI Part #MIF062), Door Bellows (RPI Part #TUB019), Rubber Boot For Door Switch (RPI Part #TUB072), Door Gasket (RPI Part #TUG021), Air Jet Valve (Red Top) (RPI Part #TUJ033), and a PM Check List

**Models:** 2340EK/MK




**RPI Part #TUK122**  
**OEM Part #02610020**  
(Door Gasket Only)

### STERILIZER PM KIT

- Includes Fill/Vent Mesh Chamber Filter (RPI Part #MIF062), Door Bellows (RPI Part #TUB019), Rubber Boot For Door Switch (RPI Part #TUB072), Door Gasket (RPI Part #TUG001), and Air Jet Valve (Red Top) (RPI Part #TUJ033), and a PM Check List

**Models:** 1730EK/MK/MKV




**RPI Part #TUK128**  
**OEM Part #02610023**  
(Door Gasket Only)

### STERILIZER PM KIT

- Includes Fill/Vent Mesh Chamber Filter (RPI Part #MIF062), Door Bellows (RPI Part #TUB019), Rubber Boot For Door Switch (RPI Part #TUB072), Door Gasket (RPI Part #TUG003), Air Jet Valve (Black Top) (RPI Part #TUJ034), and a PM Check List

**Models:** 2540E/M




**RPI Part #TUK123**  
**OEM Part #02610005**  
(Door Gasket Only)

### STERILIZER PM KIT

- Includes Fill/Vent Mesh Chamber Filter (RPI Part #MIF062), Door Bellows (RPI Part #TUB019), Rubber Boot For Door Switch (RPI Part #TUB072), Door Gasket (RPI Part #TUG002), Air Jet Valve (Black Top) (RPI Part #TUJ034), and a PM Check List

**Model:** 2340M (Serial #8805 and below).




**RPI Part #TUK130**  
**OEM Part #02610023**  
(Door Gasket Only)

### STERILIZER PM KIT

- Includes Fill/Vent Mesh Chamber Filter (RPI Part #MIF062), Door Bellows (RPI Part #TUB019), Rubber Boot For Door Switch (RPI Part #TUB072), Door Gasket (RPI Part #TUG003), Air Jet Valve (Red Top) (RPI Part #TUJ033), and a PM Check List

**Models:** 2540EK/MK




**RPI Part #TUK124**  
**OEM Part #02610118**  
(Door Gasket Only)

### STERILIZER PM KIT

- Includes Fill/Vent Mesh Chamber Filter (RPI Part #MIF062), Door Bellows (RPI Part #TUB019), Rubber Boot For Door Switch (RPI Part #TUB072), Door Gasket (RPI Part #TUG021), Air Jet Valve (Black Top) (RPI Part #TUJ034), and a PM Check List

**Models:** 2340E & 2340M (Serial #8806 and above).




**RPI Part #TUK131**  
**OEM Part #02610019**  
(Door Gasket Only)

### STERILIZER PM KIT

- Includes Fill/Vent Mesh Chamber Filter (RPI Part #MIF062), Door Bellows (RPI Part #TUB019), Rubber Boot For Door Switch (RPI Part #TUB072), Door Gasket (RPI Part #TUG074), Air Jet Valve (Black Top) (RPI Part #TUJ034), and a PM Check List

**Models:** 3870E/M

## PM Kits to fit Tuttnauer EZ9, EZ10, EZ10k, 1730, 2340, 2540 & 3870 series - continued




**RPI Part #TUK125**  
**OEM Part #02610118**  
 (Door Gasket Only)

### STERILIZER PM KIT

- Includes Fill/Vent Mesh Chamber Filter (RPI Part #MIF062), Door Bellows (RPI Part #TUB019), Rubber Boot For Door Switch (RPI Part #TUB072), Door Gasket (RPI Part #TUG021), and Air Jet Valve (Black Top) (RPI Part #TUJ034), HEPA Filter (RPI Part #RPF644), and a PM Check List

**Models:** EZ9 & 2340EA




**RPI Part #TUK129**  
**OEM Part #02610023**  
 (Door Gasket Only)

### STERILIZER PM KIT

- Includes Fill/Vent Mesh Chamber Filter (RPI Part #MIF062), Door Bellows (RPI Part #TUB019), Rubber Boot For Door Switch (RPI Part #TUB072), Door Gasket (RPI Part #TUG003), Air Jet Valve (Red Top) (RPI Part #TUJ033), HEPA Filter (RPI Part #RPF644), and a PM Check List

**Models:** EZ10, EZ10k, 2540EA/EKA/EHS




**RPI Part #TUK127**  
**OEM Part #02610118**  
 (Door Gasket Only)

### STERILIZER PM KIT

- Includes Fill/Vent Mesh Chamber Filter (RPI Part #MIF062), Door Bellows (RPI Part #TUB019), Rubber Boot For Door Switch (RPI Part #TUB072), Door Gasket (RPI Part #TUG021), Air Jet Valve (Red Top) (RPI Part #TUJ033), HEPA Filter (RPI Part #RPF644), and a PM Check List

**Model:** 2340EKA



**RPI Part #TUK132**  
**OEM Part #02610019**  
 (Door Gasket Only)


### STERILIZER PM KIT

- Includes Fill/Vent Mesh Chamber Filter (RPI Part #MIF062), Door Bellows (RPI Part #TUB019), Rubber Boot For Door Switch (RPI Part #TUB072), Door Gasket (RPI Part #TUG074), Air Jet Valve (Black Top) (RPI Part #TUJ034), HEPA Filter (RPI Part #RPF644), and a PM Check List

**Models:** 3870EA/EHS

## PM Kits to fit Pelton & Crane Delta AE8, XL8, Q8, AF10, XL10 & Q10

**RPI Part #PCK231**  
**OEM Part #1539241**  
 (Door Gasket Only)




### STERILIZER PM KIT

- Includes Compressor Filter (RPI Part #SCF002), Door Gasket (RPI Part #PCG065), Filter (Fill Line) (RPI Part #PCF009), Filter In-Line Assembly (RPI Part #PCF233), HEPA Filter (RPI Part #RPF644), Reservoir Filter Kit (RPI Part #PCK167), Spacer (Door) (RPI Part #PCS067), and a PM Check List

**Models:** AE8, XL8 & Q8

**RPI Part #PCK232**  
**OEM Part #1539407**  
 (Door Gasket Only)



### STERILIZER PM KIT


- Includes Compressor Filter (RPI Part #SCF002), Door Gasket (RPI Part #PCG066), Filter (Fill Line) (RPI Part #PCF009), Filter In-Line Assembly (RPI Part #PCF233), HEPA Filter (RPI Part #RPF644), Reservoir Filter Kit (RPI Part #PCK167), Spacer (Door) (RPI Part #PCS067), and a PM Check List

**Models:** AF10, XL10 & Q10

More PM Kits to fit Pelton & Crane on page 4!

# Table Top Sterilizers


## PM Kits to fit Pelton & Crane Magnaclave, OCM, OCR, OCR+, Sentry & Validator



**RPI Part #PCK216**  
**OEM Part #004497**  
 (Door Gasket Only)

### STERILIZER PM KIT


- Includes Door Gasket (RPI Part #PCG019), Air Valve O-ring (RPI Part #RPO354), Bellows (RPI Part #PCB001), Filter with Cotter Pin (RPI Part #PCK226), Filter (Fill Line) (RPI Part #PCF009), and a PM Check List
- Model:** Magnaclave



**RPI Part #PCK220**  
**OEM Part #1539241**  
 (Door Gasket Only)

### STERILIZER PM KIT


- Includes Door Gasket (RPI Part #PCG065), Filter (Fill Line) (RPI Part #PCF009), Air Valve O-ring (RPI Part #RPO354), Bellows (RPI Part #PCB001), Spacer (Door) (RPI Part #PCS067), and a PM Check List
- Model:** Validator 8



**RPI Part #PCK217**  
**OEM Part #001187**  
 (Door Gasket Only)

### STERILIZER PM KIT


- Includes Door Gasket (RPI Part #PCG020), Filter (Fill Line) (RPI Part #PCF009), Air Valve O-ring (RPI Part #RPO354), Bellows (RPI Part #PCB001), and a PM Check List
- Model:** Sentry



**RPI Part #PCK221**  
**OEM Part #1539407**  
 (Door Gasket Only)

### STERILIZER PM KIT


- Includes Door Gasket (RPI Part #PCG066), Filter (Fill Line) (RPI Part #PCF009), Air Valve O-ring (RPI Part #RPO354), Bellows (RPI Part #PCB001), Spacer (Door) (RPI Part #PCS067), and a PM Check List
- Model:** Validator 10



**RPI Part #PCK218**  
**OEM Part #004014**  
 (Door Gasket Only)

### STERILIZER PM KIT


- Includes Door Gasket (RPI Part #PCG014), Filter (Fill Line) (RPI Part #PCF009), Air Valve O-ring (RPI Part #RPO354), Bellows (RPI Part #PCB001), and a PM Check List
- Model:** OCM



**RPI Part #PCK222**  
**OEM Part #1539241**  
 (Door Gasket Only)

### STERILIZER PM KIT


- Includes Door Gasket (RPI Part #PCG065), Filter (Fill Line) (RPI Part #PCF009), Spacer (Door) (RPI Part #PCS067), Reservoir Filter Kit (RPI Part #PCK167), and a PM Check List
- Models:** Validator Plus 8



**RPI Part #PCK219**  
**OEM Part #004341**  
 (Door Gasket Only)

### STERILIZER PM KIT

- Includes Door Gasket (RPI Part #PCG015), Filter (Fill Line) (RPI Part #PCF009), Air Valve O-ring (RPI Part #RPO354), Bellows (RPI Part #PCB001), and a PM Check List
- Model:** OCR & OCR+



**RPI Part #PCK223**  
**OEM Part #1539407**  
 (Door Gasket Only)

### STERILIZER PM KIT

- Includes Door Gasket (RPI Part #PCG066), Filter (Fill Line) (RPI Part #PCF009), Spacer (Door) (RPI Part #PCS067), Reservoir Filter Kit (RPI Part #PCK167), and a PM Check List
- Models:** Validator Plus 10