


















New Parts to fit Midmark® • Ritter® Sterilizers

M7 & 7 SpeedClave®

MIDMARK® • RITTER®

  <p>RPI Part #RCV112 OEM Part #014-0236-00</p> <p>SOLENOID VALVE (FILL/VENT) Fits: Mounts on Chassis Models: M7-020, -022 All Serial Numbers with V Prefix</p>	  <p>RPI Part #RCK113 OEM Part #014-0236-11</p> <p>REPAIR KIT (FILL/VENT) Fits: Solenoid Valve (Fill/Vent) Models: M7-020, -022 All Serial Numbers with V Prefix</p>	  <p>RPI Part #RCC114 OEM Part #014-0236-10</p> <p>COIL (FILL/VENT) Fits: Solenoid Valve (Fill/Vent) Models: M7-020, -022 All Serial Numbers with V Prefix</p>
<p>Need Repair Kit? →</p>		<p>Need Coil? →</p>

  <p>RPI Part #RCS115 OEM Part #015-1791-00</p> <p>FILL/VENT SWITCH Fits: Front Panel Models: M7-020 thru -022 All Serial Numbers with V Prefix</p>	 <p>RPI Part #RPT501 OEM Part # (No OEM # Available)</p> <p>DOUBLE ENDED SPANNER WRENCH Fits: Solenoid Valve (Fill/Vent) Models: M7-020 thru -022 All Serial Numbers with V Prefix</p>	 <p>RPI Part #ILS014 OEM Part # (No OEM # Available)</p> <p>SPANNER NUT Fits: Solenoid Valve (Fill/Vent) Models: M7-020 thru -022 All Serial Numbers with V Prefix</p>
---	--	--

 <p>RPI Part #MIF062 New OEM Part #002-0360-00 Old OEM Part #016-0790-00</p> <p>FILL/VENT MESH CHAMBER FILTER (Stainless Steel) Fits: Chamber Models: M7-020 thru -022 All Serial Numbers with V Prefix</p>	  <p>RPI Part #RPH674 OEM Part #015-1259-04</p> <p>FUSE HOLDER Fits: Rear Panel Models: M7-020, -022 All Serial Numbers with V Prefix</p>	  <p>RPI Part #RPF679 OEM Part #015-0346-35</p> <p>FUSE (12A, 250VAC) - FAST ACTING (5/pkg) Fits: Fuse Holder (RPI Part #RPH674) Models: M7-020, -022 All Serial Numbers with V Prefix</p>
--	--	---










<table border="1" style="width: 100%;"> <tr> <td style="width: 50%; text-align: center;"> <p>Plug Type</p>  </td> <td style="width: 50%; text-align: center;"> <p>Connector Type</p>  </td> </tr> </table> <p>RPI Part #PCC003 OEM Part #H98159</p> <p>INDUSTRIAL GRADE POWER CORD Fits: Lower Chassis Models: M7 & 7 SpeedClave® All Models</p>	<p>Plug Type</p> 	<p>Connector Type</p> 	 <p>RPI Part #RPB167 OEM Part #H98160</p> <p>STRAIN RELIEF BUSHING (6/pkg) Fits: Power Cord on Lower Chassis Models: M7 & 7 SpeedClave® All Models</p>	 <p>RPI Part #RPT172 OEM Part # (No OEM # Available)</p> <p>STRAIN RELIEF BUSHING PLIERS Fits: Strain Relief Bushing Models: M7 & 7 SpeedClave® All Models</p>
<p>Plug Type</p> 	<p>Connector Type</p> 			

Table Top Sterilizers

M7 & 7 SpeedClave®

RPI Part #RCB088
OEM Part #H96772

T-BOLT
Fits: Door Assembly
Models: M7 & 7 SpeedClave®
All Models



RPI Part #RCG104
OEM Part #053-1579-00

FIBER GASKET (5/pkg)
Fits: Door Assembly
Models: M7 & 7 SpeedClave®
All Models



RPI Part #RCW009
OEM Part #H96771 & 3059

DOOR WASHER
Fits: T-Bolt on Door Assembly
Models: M7 & 7 SpeedClave®
All Models




RPI Part #RCN094
OEM Part #H96773

CAP NUT
Fits: T-Bolt in Door Assembly
Models: M7 & 7 SpeedClave®
All Models



RPI Part #RPA032
OEM Part #042-0024-02

THREADLOCKER 242 (5/pkg)
Fits: Cap Nut and T-Bolt on the Door Assembly
Models: M7




RPI Part #RCS108
OEM Part #H106330

DOOR SPRING (5/pkg)
Fits: Door Assembly
Models: M7 & 7 SpeedClave®
All Models




RPI Part #RPK431
OEM Part #002-0653-00

DRAIN TUBE ASSEMBLY KIT
Fits: Front Panel
Models: M7 (-011 thru -016) & (-020 thru -022)



RPI Part #RPF427
OEM Part #053-0990-01

MALE QUICK CONNECT FITTING (GREY)
Fits: Drain Tube Assembly Kit (RPI Part #RPK431)
Models: M7 (-011 thru -016) & (-020 thru -022)



RPI Part #RPO304
OEM Part # (No OEM # Available)
Industry Part #AS568-011


O-RING (12/pkg)
Fits: Male Quick Connect Fitting (Grey) (RPI Part #RPF427)
Models: M7 (-011 thru -016) & (-020 thru -022)



SHOWN
ACTUAL
SIZE

RPI Part #RPF426
OEM Part #053-0990-00

FEMALE QUICK CONNECT FITTING (GREY)
Fits: Front Panel
Models: M7 (-011 thru -016) & (-020 thru -022)




RPI Part #RPH673
OEM Part #040-0010-125

SELF DRILLING SCREW (6/pkg)
Fits: Cover
Models: M7 & 7 SpeedClave®
All Models



RPI Part #RPH285
OEM Part #053-0391-03 & 053-0391-02

DRAIN HOSE
Fits: Fits: Male Quick Connect Fitting (Grey) (RPI Part #RPF427)
Models: M7 (-011 thru -016) & (-020 thru -022)



Sold by the foot; Reel not included

RPI Part #RPC286
OEM Part # (No OEM # Available)

KWIK™ CLAMP (6/pkg)
Fits: Drain Hose/Inner Tubing
Models: M7 (-011 thru -016) & (-020 thru -022)



SERVICE TIP


The efficiency of a sterilizer itself should be monitored periodically to determine if proper sterilization is taking place. And one of the ways to do this is to check that the sterilizer is performing at the proper heat temperature.

A good tool to check the temperature inside a sterilizer is the RPI Max Register Thermometer. This thermometer registers the highest temperature it has been exposed to and holds that temperature until reset. To reset, simply shake it until the mercury inside reaches ambient temperature.

For more ways to monitor a sterilizer for optimal performance, please see the "Tech Talk" article in our September, 2008 newsletter.

RPI Part #RPT113
OEM Part # (No OEM # Available)

MAX REGISTER THERMOMETER
• Use to test temperature inside sterilizer



Look for more parts to fit M7 & 7 SpeedClave Sterilizers coming soon!