









DIN/UNSHIELDED LEADWIRE SETS FOR VETERINARY USE

RPI PART #	DESCRIPTION	# OF LEADS	LENGTH	EQUIPMENT END	PATIENT END
LVM701	3 LEADWIRE SET - 24"DN/AL • Reusable	3	24"	DIN/UNSHIELDED .060 SAFETY SOCKET 	ALLIGATOR CLIP  WHITE RED BLACK
LVM702	3 LEADWIRE SET - 24"DN/AL • Reusable	3	24"	DIN/UNSHIELDED .060 SAFETY SOCKET 	ALLIGATOR CLIP  WHITE GREEN BLACK
LVM703	5 LEADWIRE SET - 24"DN/AL • Reusable	5	24"	DIN/UNSHIELDED .060 SAFETY SOCKET 	ALLIGATOR CLIP  RED WHITE GREEN BLACK BROWN
LVM704	7 LEADWIRE SET - 24"DN/AL • Reusable	7	24"	DIN/UNSHIELDED .060 SAFETY SOCKET 	ALLIGATOR CLIP  ORANGE BLUE BROWN BLACK GREEN WHITE RED