



RPI PART #	DESCRIPTION	AIR TECHNIQUES	APOLLO/ MIDMARK	DCI/MDT McKESSON	DENTALEZ®/ CUSTOM AIR/ RAMVAC	JUN-AIR	MATRIX/ MIDMARK	MIDMARK	TECH WEST/ TURBINE,IND.
AIR INTAKE FILTERS									
CME004	REPLACEMENT FILTER ELEMENT (FITS RPI PART #CMK003 ONLY)		PFM50955						
CMF002	AIR INTAKE FILTER (METAL STYLE)	80050	PFM50955						
CMF011	FOAM ELEMENT	60043-1							
CMF013	FOAM ELEMENT						77007079		AAF-100
CMF014	FOAM ELEMENT			3-08-0632-10					
CMF073	AIR INTAKE FILTER					JA5413210			
CMF075	INTAKE FILTER					JA5412400			
CMK003	FRESH AIR INTAKE FILTER KIT Includes Replaceable Filter Element (RPI Part #CME004)		PFM50955					ACP90505	
COOLING FANS									
CMF019	COOLING FAN (115VAC)	41060	SEA95535	3-08-0578-10	5451400	JA5445400			CF-115
CMF020	COOLING FAN (230VAC)	80516	SEA95540	3-08-0563-10	5451500	JA5445000			CF-230
COMPRESSOR OIL									
CML001	COMPRESSOR OIL (1) 10 oz. Bottle	60051	ACA85921	3-08-0084-30	64651001		77000220		CO-100
CML045	COMPRESSOR OIL (12) 16 oz. Bottles	60051	ACA85922	3-08-0084-30	64651001		77000220		CO-100
DESICCANT									
CMD091	DESICCANT	60010	RAM80105						
CMD093	DESICCANT		RAM80100						MDS-200
CMD094	DESICCANT								MDS-100
FILTER ASSEMBLIES & COMPONENTS									
CMA021	COALESCING FILTER ASSY (5 OZ.) Wilkerson brand M16		ACA70305				77005010	62985200	
CMA022	COALESCING FILTER ASSY (10 OZ.) Wilkerson brand M26		ACA85215		645-45010	JA4081500			CFA-375
CMA068	AIR FILTER ASSEMBLY								TAF-100
CMA070	AIR FILTER ASSEMBLY								TAF-200
CMA072	AIR FILTER ASSEMBLY				004006				
CMA082	PARTICLE FILTER ASSEMBLY Norgren brand F73	85617							
CMA084	PARTICLE FILTER ASSEMBLY Wilkerson brand F18	86197							
CMB087	MOUNTING BRACKET Wilkerson brand F16 & M16 style		No OEM Part # Available				No OEM Part # Available	No OEM Part # Available	
CMB088	MOUNTING BRACKET Wilkerson brand F26 & M26 style		No OEM Part # Available				No OEM Part # Available	No OEM Part # Available	
CMB099	MOUNTING BRACKET Norgren brand F73 style	No OEM Part # Available							
CMB100	MOUNTING BRACKET Wilkerson brand F18 & M18 style	No OEM Part # Available							
CMC086	BLANK COVER Wilkerson brand M16 & M26 style		No OEM Part # Available				No OEM Part # Available	No OEM Part # Available	
CME043	COALESCING FILTER ELEMENT Wilkerson brand (Type B)				004016				
CME044	COALESCING FILTER ELEMENT Wilkerson brand (Type C)				003846				
CME067	AIR FILTER ELEMENT				003965				TIF-100
CME069	AIR FILTER ELEMENT								TIF-200
CME071	AIR FILTER ELEMENT						77001541		
CME080	AIR FILTER ELEMENT				003965				
CME083	PARTICLE FILTER ELEMENT Norgren brand F73 style	85602							
CME085	PARTICLE FILTER ELEMENT Wilkerson brand FF18 & M18 style	86193							
CME089	CARBON FILTER ELEMENT Wilkerson brand (Type D)		ACA70305*				77005010*	62985200*	
CME090	CARBON FILTER ELEMENT Wilkerson brand (Type D)		ACA85215*		645-45010*	JA4081500*			CFA-375
CMF016	AUTO FLOAT DRAIN Wilkerson brand		ACA85210		645-45010				
CMIO15	DIFFERENTIAL PRESSURE INDICATOR Wilkerson brand		AAP85210						CI-100



RPI PART #	DESCRIPTION	AIR TECHNIQUES	APOLLO/ MIDMARK	DCI/MDT McKESSON	DENTALEZ [®] / CUSTOM AIR/ RAMVAC	JUN-AIR	MATRIX/ MIDMARK	MIDMARK	TECH WEST/ TURBINE,IND.
FILTER ASSEMBLIES & COMPONENTS (continued)									
CMIO95	DIFFERENTIAL PRESSURE INDICATOR Wilkerson brand B18, F18 & M18 style	86192							
CMIO96	DIFFERENTIAL PRESSURE INDICATOR Norgren brand F73 style	86192							
CMK005	COALESCING FILTER ELEMENT Wilkerson brand (Type B)		ACA86050		648-00209	648-00209	77005010		
CMK006	COALESCING FILTER ELEMENT Wilkerson brand (Type B)		ACA86055		645-45010	JA40815000			CFE-375
CMK024	COALESCING FILTER ELEMENT Wilkerson brand (Type C)		ACA86050		648-00209	648-00209	77005005		
CMK025	COALESCING FILTER ELEMENT Wilkerson brand (Type C)		ACA86055		645-45010	JA40815000			CFE-375
CMK026	COALESCING FILTER ELEMENT Generic brand (Type C)		ACA86050		645-00209	648-00209	77005005		
CMK027	COALESCING FILTER ELEMENT Generic brand (Type C)		ACA86055		645-45010	JA40815000			CFE-375
CMK076	COMPRESSOR PM KIT				004016				
RPF651	ADAPTER FITTING Fits Auto Float Drain		ACA85210		645-45010				
MOISTURE INDICATORS									
CMIO77	MOISTURE INDICATOR	60013	ACA86965						
CMIO78	MOISTURE INDICATOR				004021				MI-100
CMIO79	MOISTURE INDICATOR						77005140		
CMIO81	MOISTURE INDICATOR			3-08-1019-20					
PRESSURE GAUGES									
CMG048	PRESSURE GAUGE (160 PSI)	60011	PCA70400				77005012		CPG-250
CMG049	PRESSURE GAUGE (300 PSI) 1/4 MPT bottom mount	No OEM Part # Available	No OEM Part # Available	No OEM Part # Available	No OEM Part # Available	No OEM Part # Available	No OEM Part # Available	No OEM Part # Available	No OEM Part # Available
CMG050	PRESSURE GAUGE (230 PSI)					JA5130000			
PRESSURE SWITCHES									
CMS051	PRESSURE SWITCH		ECS10445						
CMS052	PRESSURE SWITCH		ECS10455						
CMS053	PRESSURE SWITCH								PCS-3
CMS054	PRESSURE SWITCH					JA5074800			
CMS055	PRESSURE SWITCH					JA5075200			
CMS056	PRESSURE SWITCH	85465							
CMS057	PRESSURE SWITCH	60016							
RPO303	O-RING (12/pkg)					No OEM Part # Available			
RPT670	TUBING (6mm OD Black Poly.) (Sold by the foot)					No OEM Part # Available			
PURGE FILTERS & MUFFLERS									
CMF012	FOAM INSERT/PURGE FILTER								FA-100
CMM017	PURGE MUFFLER	80330							AS-250
CMM023	PURGE MUFFLER								PM-500
PRESSURE RELIEF VALVES, SAFETY VALVES & UNLOADER VALVE									
CMV059	PRESSURE RELIEF VALVE	No OEM Part # Available	PVV05526	No OEM Part # Available					
CMV060	PRESSURE RELIEF VALVE		No OEM Part # Available	No OEM Part # Available					No OEM Part # Available
CMV046	SAFETY VALVE (125 PSI)	60017	PW50560	3-08-0813-00			77005020		POVA-100
CMV047	SAFETY VALVE (135 PSI)					JA5417000			
CMV058	UNLOADER VALVE					JA6973300			
RELAY (UNIVERSAL)									
CMR032	RELAY (UNIVERSAL)	60501 62501 64501 66501	HCA30652 HCA30653 HCA30654 HCA30655		645011010* 645011011*		77001344* 77001346*	77001344* 77001346*	PR-115 PR-208
RUN CAPACITORS									
CMC031	CAPACITOR (20µf, 440V, 50/60Hz.)	66516	HCA30622				77001346*	77001346*	RCC-208
CMC034	CAPACITOR (10µf, 370V, 50/60Hz.)				645011011*				
CMC037	CAPACITOR (40µf, 370V, 50/60Hz.)	62503	HCA30620		645011010*		77001346*	77001346*	RCC-115



RPI PART #	DESCRIPTION	AIR TECHNIQUES	APOLLO/ MIDMARK	DCI/MDT McKESSON	DENTAL EZ [®] / CUSTOM AIR/ RAMVAC	JUN-AIR	MATRIX/ MIDMARK	MIDMARK	TECH WEST/ TURBINE,IND.
SIGHT GLASSES & SIGHT GLASS TOOL									
CMK097	SIGHT GLASS	60100	ACA85322				77000812		OSG-100
CMK098	SIGHT GLASS KIT	60100	ACA85322				77000812		OSG-100
RPT705	SIGHT GLASS TOOL	No OEM Part # Available	ACA00050				No OEM Part # Available	No OEM Part # Available	No OEM Part # Available
START BOXES									
CMA029	START BOX (1HP, 208/230V)	66515	SPD95291				77001346	77001346	CRA-208
CMA033	START BOX (1HP, 230V)	66500	SPD95290		645011011		77001344		CRA-230
CMA035	START BOX (1HP, 115V)	62500	SPD95280		645011010		77001345		CRA-115
CMA039	START BOX (3/4HP, 208/230V)	64515	SPD95285 SPD95286						
CMA041	START BOX (3/4HP, 115V)	60500	SPD95275						
START CAPACITORS									
CMC030	CAPACITOR (145-174µf, 250V, 50/60Hz.)	64502	HCA30625				77001344* 77001345* 77001346*	77001344* 77001345* 77001346*	SCC-230
CMC036	CAPACITOR (430-516µf, 115V, 50/60Hz.)	62502	HCA30628				77001345*	77001345*	SCC-115
CMC042	CAPACITOR (540-648µf, 115V, 50/60Hz.)	60502	HCA30626						
THERMAL PROTECTOR ASSEMBLIES									
CMP061	THERMAL PROTECTOR ASSEMBLY	66504	HCA30615		No OEM Part # Available		61894802	61894802	KOPC-230
CMP062	THERMAL PROTECTOR ASSEMBLY	62504	HCA30610		No OEM Part # Available		61894801	61894801	KOPC-230
CMP063	THERMAL PROTECTOR ASSEMBLY	64503	HCA30605						
CMP064	THERMAL PROTECTOR ASSEMBLY	64503	HCA30600						

* The OEM does not list an individual part number, however the part is included in the OEM Part Number noted.